# Ecosite Bond

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## English

### Product description

Ecosite Bond is a light-curing single-component bonding agent. The material can be used as an adhesive on enamel and dentine with the self-etch, total-etch or selective-etch technique in combination with light-curing filling materials.

### Indications

Bonding agent for direct restorations of all cavity classes with light-curing composites, compomers, and ormocers.

### Contra-indications

- Do not use the material where there are known allergies to any of the components or contact allergies.
- Do not use directly on the exposed pulp.
- Do not use the material if a relatively dry working area or the recommended application technique is not possible.

### Interactions

Moisture, materials containing eugenol and oily air may inhibit polymerization at the contact area.

### Side effects

There are no known side effects to date. Hypersensitivity to components of the material cannot be ruled out in individual cases.

### Basic safety information

- For dental use only!
- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid eye contact! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Avoid contamination of the bottle content from unclean disposable brushes.

### Notes for use

- ▶ Light-curing units should have an output of 450 nm and should be checked regularly. The light intensity should be a minimum of 400 mW/cm<sup>2</sup>. Place the light as close as possible to the material. Observe the manufacturer's instructions.
- ▶ In order to prevent premature polymerization due to ambient light, the material must be added to the mixing palette immediately prior to application.
- ▶ Ensure that the etched or bonded surface does not become contaminated prior to each subsequent working step.
- ▶ Do not remove the oxygen inhibition layer caused by the application because this is important for the bonding process with the material applied on top.
- ▶ Observe the manufacturer's instructions for all materials that are used with Ecosite Bond.
- ▶ Do not subject the dentine to excessive drying, particularly after etching.

### Recommended use

#### Prepare tooth structure and etch if necessary

1. Prepare the cavity in accordance with the general rules for adhesive technique. The use of a rubber dam is recommended. In deep cavities, protect areas close to pulp with a suitable liner.
2. Where necessary, etch as follows:
  - ▶ Selective-etch technique: Etch affected enamel with 37% phosphoric acid gel. The response time of the etching gel on the enamel should be 20 to 60 s. The etched enamel should have a chalky white appearance. If this is not the case, etching must be repeated. Accidental etching of dentine has no effect on the bonding values, but the response time of the etching gel on dentine should not exceed 15 s.
  - ▶ Total-etch technique: Etch dentine and enamel with 37% phosphoric acid gel. Start the etching process with the enamel and then expand this to include the dentine. The response time of the etching gel on the enamel should be 20 to 60 s. The etched enamel should have a chalky white appearance. If this is not the case, etching must be repeated. The response time on the dentine must not exceed 15 seconds.
3. Rinse with water for at least 15 seconds. Remove large areas of excess water.

4. After drying, an even, moist layer must remain on the etched dentine surface requiring treatment. If necessary, wet the dentine surface using a lightly moistened micro brush.

#### Apply adhesive to the tooth structure

1. Place 2–3 drops of Ecosite Bond on a mixing palette. Using a micro brush, rub the bond in for 10 s before leaving to work for 20 s.
2. Blow the materials for 5 s to remove the solvent.

**Note:** Ensure that enough material remains on the tooth.

3. Light-cure the bond with an LED polymerization device or a halogen polymerization device for 10 s.
4. Apply the filling material to the cured Ecosite Bond.

#### Composition

Dental resins, ethanol, water, additives and catalysts

#### Storage and shelf life

- Store in a dry place at room temperature (2 to 25 °C/36 to 77 °F)!
- Close immediately after use so that the material does not cure in the ambient light.
- Do not use after the expiration date.

**Note:** A glossary of all symbols displayed is at: <http://www.dmg-america.com/literature>

**Caution:** US Federal Law restricts this device to sale by, or on the order of a dentist, or other practitioner licensed by law of the state in which he or she practices to use or order the use of this device.

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