# VITIQUE

esthetic resin cement

**DMG**

DE – Gebrauchsinformation	4
EN – Instructions for use	11
FR – Mode d'emploi	18
IT – Istruzioni per l'uso	25
ES – Instrucciones de uso	33
PT – Instruções de utilização	40
NL – Gebruiksaanwijzing	47
DA – Brugsanvisning	55
SV – Bruksanvisning	62
PL – Instrukcja użycia	69
RU – Информация по применению	76
TR – Kullanım bilgileri	84
ZU – 使用信息	91
KO – 사용 정보	97
JA – 使用説明書	104

## **Produktbeschreibung**

---

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## English

### Product Description

The Vitique Veneer Cementation System is used for the permanent, highly esthetic cementation of veneers made from ceramic or composite. The Vitique Veneer Cementation System contains the following components:

#### I. Vitique Veneer Cementation Kit

**Vitique Try-In:**

Using the Vitique Try-In pastes, the final shade in the mouth can be accurately simulated and the correct shade of Vitique Esthetic Resin Cement can be determined.

**Vitique Esthetic Resin Cement:**

Vitique Esthetic Resin Cement is a highly esthetic composite luting cement. Vitique Esthetic Resin Cement is light-curing and X-ray opaque.

#### II. Vitique Accessory Kit:

**Vitique Silane:**

Vitique Silane is a bonding agent between silicate ceramics and composite luting cements.

**DMG Etching Gel:**

DMG Etching Gel is a universal etching gel (37% phosphoric acid gel) for enamel and dentine etching.

**TECO:**

TECO is a light-curing bonding agent for enamel and dentine for use with Vitique Esthetic Resin Cement and other light-curing composites.

## Indications

---

### **Vitique Try-In**

Paste for shade simulation.

### **Vitique Esthetic Resin Cement**

Highly esthetic cementation of silicate ceramic and composite restorations with low film thickness (< 2.0 mm) that can only be cemented through light-curing due to their high transparency (inlays, onlays and veneers).

### **Vitique Silane**

Bonding agent between etched silicate ceramics and composite luting cements.

### **DMG Etching Gel**

Enamel and dentine etching.

### **TECO**

Light-curing enamel and dentine bonding agent in conjunction with the Total Etch technique for use with light-curing composites (e.g. Vitique Esthetic Resin Cement) for indirect restorations (e.g. veneers) and direct restorations.

## Contraindications

---

- ▶ Do not use Vitique Veneer Cementation System in the case of known allergies to one of the ingredients.
- ▶ Do not apply Vitique Silane to exposed dentine or enamel.
- ▶ Use TECO only after etching the tooth substance beforehand.
- ▶ Do not use TECO and Vitique Esthetic Resin Cement directly on exposed pulp. If necessary, use suitable pulp protection.
- ▶ Do not use Vitique Esthetic Resin Cement if a dry working area on the prepared tooth is not possible. DMG recommends the use of a rubber dam.

## Basic safety instructions

---

- ▶ For dental use only!
- ▶ Keep out of reach of children!
- ▶ Avoid contact with the skin and mucous membranes. In case of accidental contact, wash the area thoroughly with water and soap and get medical attention if necessary.
- ▶ Avoid eye contact. In case of accidental eye, wash eyes with large amounts of water and get medical attention. Risk of burns!

- ▶ DMG Etching Gel can cause chemical burns. In the event of a chemical burn, seek advice from a doctor.

## Side effects

---

There are no known side effects to date. Isolated cases of hypersensitivity to components of the materials cannot be excluded.

## Interactions with other materials

---

Avoid products that contain eugenol, moisture, and air that contains oil since they can impair the curing of Vitique Esthetic Resin Cement and TECO.

## Time line

---

Mixing time Vitique Silane	15 sec
Response time Vitique Silane	10 sec
Response time DMG Etching Gel	Enamel: 20 – 60 sec Dentine: max. 15 sec
Work-in time TECO	20 sec
Light-curing TECO	10 sec
Light-curing for removal of excess Vitique Esthetic Resin Cement	1 – 3 sec
Light-curing Vitique Esthetic Resin Cement	40 sec. per side

## Instructions for use


---

- ▶ Light-curing units should have an output of 450 nm and should be checked regularly. The light intensity should be a minimum of 400 mW/cm<sup>2</sup>. Place the light as close as possible to the material. Observe the manufacturer's instructions for the light unit.
- ▶ To prevent premature curing of the materials, close all packaging immediately after use.
- ▶ The inner surfaces of silicate ceramic veneer must be etched and silanized prior to cementation. The inner surfaces of composite veneer should be sandblasted prior to cementation. Follow the manufacturer's instructions.

## Recommended use

---

### Flow Chart

1. After removing the temporary prosthesis, clean the prepared core with fluoride-free polishing paste; in doing so, remove residues of the temporary cementation material, rinse, and dry with oil-free air. Try in the restoration.
2. The cement shade is selected prior to placing the rubber dam. The shade of the luting cement is determined using the Vitique Try-In paste by applying it on the fitting surface of the veneer. If correcting the selected shade of the cement, the core and preparation are to be cleaned with water and the try-in is to be repeated with a suitable shade.
3. Remove the restoration. Clean restoration and prepared tooth thoroughly with water spray.

**Note:** Any residues of the Vitique Try-In paste can impair the permanent cementation of the restoration.

**Note:** The silane layer of pre-silanized ceramics is destroyed by the try-in. Following contamination of the silanized surface, clean the surface with acetone or alcohol and re-silanize.

4. Ensure that the restoration has been pretreated by the manufacturer (e.g., etching the silicate ceramic with hydrofluoric acid).
5. Dry the etched ceramic with oil-free air. Mix 1 to 2 drops of Vitique Silane Adhesive and Vitique Silane Activator in a 1:1 ratio for 15 sec. Immediately

apply a thin layer and leave on for 10 sec. Then carefully blow with a light stream of air.

**Note:** After contamination with Vitique Try-In Paste, clean the restoration with water spray and dry with oil- and water-free air. Clean the surface with phosphoric acid for 15 sec, rinse, dry and apply Vitique Silane .

6. The use of a rubber dam to isolate the working area is recommended. Apply DMG Etching Gel to the prepared tooth; in doing so, start in the enamel and spread to the dentine. The response time is 20 to 60 sec on the enamel and 15 sec on the dentine. Rinse the etched areas for 15 s with water and dry with oil-free air. Do not over-dry the dentine. An even, moist layer should remain on the pre-treated dentine.
7. Work the TECO bonding agent generously into the conditioned tooth substance for 20 sec using a microbrush, remove any excess present, and carefully blow using oil-free air until an evenly moist surface remains. Then light-cure for 10 sec with a suitable light unit.
8. Apply Vitique Esthetic Resin Cement with the veneer tip in the fitting surface of the restoration without leaving any bubbles. Position the restoration accurately on the prepared tooth using gentle pressure.

**Note:** By using a tacking tip (light device with max. 2 mm diameter) which is set up in the center of the veneer and activated for 1 to 6 sec, the restoration is selectively secured in its position.

9. Remove large amounts of excess material using suitable instruments. In the case of very fragile veneer edges, the excess is carefully removed with a brush or dental gauze in the direction of the gingiva.

**Note:** Excess material can be light-cured for 1 to 3 sec using a suitable light unit and can thus be removed more easily.

10. Allow each side to harden for 40 sec. Please note that the exposure time depends on the light intensity of the polymerization device used (follow manufacturer's instructions).
11. Remove excess cured material carefully using suitable instruments and then polish the edges using suitable polishing tips.
12. Check restoration for correct occlusal contacts and adjust these, if necessary.
13. Mirror-finish polishing of the edges of the restoration can be done with a polishing paste using polishing cups or discs.

## Storage and shelf life

---

- ▶ Store in a dry place at room temperature (15 to 25 °C/59 to 77 °F) !
- ▶ Do not use after the expiration date.
- ▶ Close packaging immediately after use.

## Composition

---

### Vitique Try-In

Pyrogenic silicic acid, pigments in water-soluble polyglycol matrix

### Vitique Esthetic Resin Cement

Barium glass in a Bis-GMA based matrix from dental resins, additives, catalysts.  
Filler content: 66.6% by weight = 46 by vol. (0.02 to 3.0 µm)k

### DMG Etching Gel

o-phosphoric acid in water (37%), pyrogenic silicic acid

### Vitique Silane

Silane, water, ethanol, phosphoric acid, additives

### TECO

Bis-GMA based hydrophilic and acid resin matrix, starter, additives

## Classification

---

Vitique Esthetic Resin Cement: ISO 4049:2009 (including color resistance)

TECO: ISO 29022:2013

## Packaging

---

Vitique Veneer Cementation System	REF 213080
Vitique Veneer Cementation Kit	REF 213081
Vitique Accessory Kit	REF 213082

Vitique Esthetic Resin Cement Refill	
Shade A1	REF 213038
Shade A2.5	REF 213039

<b>Vitique Esthetic Resin Cement Refill</b>	
Shade A4	REF 213040
Shade B1	REF 213041
Shade Bleach Light	REF 213042
Shade Transparent	REF 213043
Shade White	REF 213044
Shade Pink	REF 213045

<b>Vitique Try-In Refill</b>	
Shade A1	REF 213048
Shade A2.5	REF 213049
Shade A4	REF 213050
Shade B1	REF 213051
Shade Bleach Light	REF 213052
Shade Transparent	REF 213053
Shade White	REF 213054
Shade Pink	REF 213055

<b>Vitique Silane Refill</b>	
1 bottle @ 3 ml Adhesive 1 bottle @ 3 ml Activator	REF 213118

<b>TECO Refill</b>	
1 bottle @ 5 ml TECO	REF 213321

Accessories	
100 Veneer-Tips	REF 213057

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**DMG** Chemisch-Pharmazeutische Fabrik GmbH