

00:01:00

ONE MINUTE PRODUCT MEETING

O-Bite

O-Bite is the ideal partner for perfectly fitting restorations. With a working time of 30 seconds and a setting time in the mouth of 90 seconds, this superior bite registration material offers exceptional accuracy, predictability and speed.

Specifically developed by DMG to yield more stable and precise bite registrations, its proprietary vinyl polysiloxane formulation has a Shore A Hardness of 94, providing unrivaled stability and resistance to breakage. The material can be easily trimmed and finished with a scalpel or silicone cutters

With O-Bite, there's no need for spot grinding of the occlusal height, saving considerable time and providing better patient comfort. Plus, the fresh, fruity orange scent is an added bonus for both patients and doctors.



FEATURES & BENEFITS

- › Superb accuracy
- › Outstanding stability
- › Not technique sensitive
- › Hard VPS-based bite registration material
- › Superior break resistance
- › Pleasant orange scent and neutral taste
- › Superior break resistance
- › Minimal bite resistance

PROCEDURES

- › Bite registration

COMPANION PRODUCTS

- › LuxaBite
- › StatusBlue
- › Honigum Pro

ITEM # & DESCRIPTION

999765 | O-Bite

999202 | Automix Tips (Bag of 50 tips)

LOOKING FOR MORE INFO?

Visit www.dmg-america.com/products.