



Software - Hardware - Resins

DMG DentaMile

Complete, validated solution
for the digital workflow

“The 3Demax printer is extremely fast, versatile, and accurate. The “Force Feedback” technology delivers on the printing speed without compromising accuracy. The complete solution facilitates our digital production with its ease of use and the integrated and validated workflow. The small footprint and ability to stack devices adds to its versatility and convenience.”

- Edmonds Dental Prosthetics



3Demax:

Maximum reliability, minimal footprint

Our 3D printing workflow solution include the new 3Demax 3D printer, 3Dewash and 3Decure units with optimized LuxaPrint resins, accessories and cloud-based DentaMile connect software.

The patented “Force Feedback” technology of the 3Demax accelerates the printing process by up to 50%. The validated and integrated solution utilizes advanced RFID technology to ensure reliable automatic identification of materials and eliminate possible misapplications.

3DEWASH

This fully integrated unit washes 3D printed devices using a DLP process, with pre-set programs that systematically wash the remaining non-polymerized resin from the dental application.

3DECURE

This fully integrated unit cures 3D resins with 360° illumination. The pre-set programs lock in the printed device’s mechanical properties ensuring proper curing according to individual resin specifications.



Fast Print Speeds



Easy & User-Friendly



Validated Workflow

LuxaPrint resins for 3D printing

LUXAPRINT MODEL

Indicated for the fabrication of dental models that are suitable for the production of thermoformed splints.

LUXAPRINT TRAY

Indicated for the fabrication of custom impression trays.

LUXAPRINT ORTHO

Indicated for the fabrication of custom surgical guides and splints with the highest precision.

LUXAPRINT ORTHO PLUS

Indicated for the manufacturing of orthodontic and dental appliances such as mouthguards, nightguards, splints and tooth positioners.

LUXAPRINT ORTHO FLEX

A biocompatible material for the additive production of flexible splints like bleaching trays or indirect bonding trays.

LUXAPRINT CAST

Indicated for the fabrication of residue-free burn-out molds (e.g. model and denture framework castings, crowns and bridges).

LUXAPRINT GINGIVA

Ideal for creating gingival masks to help simulate the periodontal status for dental laboratory models or master models.

DENTAMILE CONNECT SOFTWARE

The intuitive DentaMile connect your cloud-based software solution that makes the entire processing and 3D production procedure easier, faster and safer. The medical device workflow production is secure and validated, errors are minimized and new applications are made accessible without any additional technical effort.

DMG America, LLC
65 Challenger Road, Suite 340, Ridgely Park, NJ 07660
tel: 800.662.6383 email: info@dmg-america.com
web: www.dmg-america.com

